

Example TCPComm Experimental File

written by John S. McCaskill and Carla Verhaelen
Ruhr University Bochum

June 2008

This is an example file for taking complex programmed control of the Omega Machine

via the interpreted language environment provided by *Mathematica*.

This file uses the mathematica package TCPComm (ver 1.2.6) written by J.S.McCaskill.

It communicates via Mathlink with an executable written in C for TCP Communication.

This in turn communicates across the network via TCP with the external command interface

for the Ng-biopro software (leading author Uwe Tangen, RUB-BioMIP) running on the experimental machine

that is connected via the Biofox stand alone interface computer with the Omega Machine.

Users can write their own command sequences (macros) and execute these in a timed interpretative environment using this software. A simple example collects data from various subimages, plotting results interactively and controlling the laser illumination synchronized with microscope images.

Set the directory to your current one containing a copy of the TCPComm executable.

The TCPComm.m mathematica package written by J.S.McCaskill also needs to be installed in your *Mathematica/Applications* folder.

```
In[1]:= SetDirectory[ "/Users/john/Library/Mathematica/Applications/TCPComm/Build/Release" ]
```

```
Out[1]= /Users/john/Library/Mathematica/Applications/TCPComm/build/Release
```

This loads the definitions from the TCPComm package (note the backquote, necessary for context name)

```
In[2]:= Needs [ "TCPComm`" ]
```

The following list of active links gives the usage information for the functions exported from the TCPComm package.

```
In[3]:= ? TCPComm` *
```

▼ TCPComm`

imbin	imextr: act	imrow: diag	imwrite	NG	NGCo: mm: une	NGPa: lette	verbo: se
imco: mbi: ne	imread	imset: ROI	MA	NGCI: ose	NGExit	status	□
imcon: vert	imROI	imshow	MAC: mds	NGC: md	NGIns: tall	Timed: Loo: p	
imcsu: ROI	imrow: col	imsize	MAPa: lette	NGC: mds	NGlink	untool	

Setup palettes once on your machine

NGPalette[]

Sess open	Sess close	Quit	Status	BackOff	BackOn
Synch	Window	Cam par	Snap	AOTF	Filter
Pump par	Pump	List pins	Electrode	Cycle	List sensors
Sensor data	Intensity	List names	Temp	XY move	XY def
ZStage	Verbose				

MAPalette[]

Timed Loop	Set ROI	I Snap	I Convert	I Display	I Dispsub
I Color	I Write	I Read	I Size	I RowCol	I Line
I Bin	I ROI	I Comb	NG Start	NG Exit	

1. Install connection to ng-biopro

```
In[9]:= MA[NGInstall["dyck"]]
```

Command completed successfully in session Mathematica_TCP

Command completed successfully in session session1

verbose mode turned off


```
Out[9]:= {LinkObject[./TCPComm, 7, 5], {Mathematica_TCP}, {session1},
  False, {Pump use desired, Graphics are up and operational, 0}}
```

2. Experimental macro session

```
In[10]:= im1 = imconvert[NG["i"]];
```

Image[1004,1002] Thu 21 Aug 2008 13:19:36

```
In[11]:= Print[imshow[im1, Red, 0, {{0, 0}, {0, 0}}]];
```



e.g. Low bleaching acquisition with measurement field

Definitions for region extraction

```
In[13]:= size[roi_] := (roi[[2, 1]] - roi[[1, 1]] + 1) (roi[[2, 2]] - roi[[1, 2]] + 1);
```

```
In[14]:= flipROI[{{x1_, y1_}, {x2_, y2_}}, {ydim_, xdim_}] :=  
Reverse[{{x2, ydim - 1 - y1}, {x1, ydim - 1 - y2}}];
```

Macro for low bleaching

```
In[27]:= laser[j_Integer] := Module[{mim},  
  NG["a", "biofox_blue", 63];  
  Pause[0.2]; im1 = imconvert[NG["i"]];  
  NG["a", "biofox_blue", 0];  
  mim = Flatten[imextract[im1[[1]], ROImess[1]]];  
  mv[1][j] = Fold[Plus, 0, mim] / size[ROImess[1]];  
  mim = Flatten[imextract[im1[[1]], ROImess[2]]];  
  mv[2][j] = Fold[Plus, 0, mim] / size[ROImess[2]];  
  PrintTemporary["Printing image"];  
  Print[GraphicsRow[{Show[imshow[im1, White, 0, imcsuROI], PlotLabel -> j],  
    ListPlot[Array[mv[1], j]], ListPlot[Array[mv[2], j]]}]]];  
]
```

Custom version of TimedLoop macro

```
In[16]:= TimedLoop1[name_String, body_, intvlsec_, niter_Integer] := Module[{i, t, t1, nwait, d},  
  For[t = SessionTime[]; i = 1, i ≤ niter, i++,  
    t1 = SessionTime[]; Print["\n", name, "[", i, "] Time = ", t1 - t, "s"];  
    body[i];  
    If[i == niter, Break[]];  
    t1 = SessionTime[];  
    nwait = t + intvlsec * i - t1;  
    If[nwait ≤ 0, Print[nwait, " is too little time for macro timing"],  
      Print["waiting..."]; Pause[nwait]];  
]
```

Take sample image

```
In[21]:= im1 = imconvert[NG["i"]][[1]];
```

```
Image[1004,1002] Thu 21 Aug 2008 13:21:23
```

First guess, ROI for a sample measurement field

```
In[22]:= roi = ROIess[1] = ROIess[2] = imsetROI[{{300, 300}, {350, 350}}]
```


```
Out[22]= {{300, 300}, {350, 350}}
```

Optimize interactively

Optimize measurement field

```
In[23]:= Dynamic[
  {LocatorPane[Dynamic[roi], Show[imshow[im1, White, 0, {{0, 0}, {0, 0}], False],
    ImageSize -> {400, 400}], Automatic, {1, 1}],
  Dynamic[
    roi]]]
```

Optimize measurement field



```
Out[23]= {
```

```
  {{423., 78.}, {517., 200.}}
```

Store vertically flipped ROI as measurement field 1

```
In[24]:= ROIess[1] = flipROI[Round[roi], Dimensions[im1]]
```

```
Out[24]= {{421, 506}, {538, 587}}
```

Repeat optimization for second measurement field. The Store vertically flipped ROI as measurement field 2


```
In[25]:= ROIess[2] = flipROI[Round[roi], Dimensions[im1]]
```

```
Out[25]= {{423, 801}, {517, 923}}
```

```
In[28]:= NG["s 0x8000 0x8000"]; TimedLoop1["laser", laser, 15, 3]; NG["s 0x8000 0"];
```

laser[1] Time = 0.000049s


Image[1004,1002] Thu 21 Aug 2008 13:25:14



waiting...

laser[2] Time = 15.006485s


Image[1004,1002] Thu 21 Aug 2008 13:25:14



waiting...

laser[3] Time = 30.009190s

Image[1004,1002] Thu 21 Aug 2008 13:25:14



3. To finish session

```
In[29]:= MA[NGExit["dyck"]]
```

verbose mode turned on

No return parameter, returned with stype è

```
Out[29]:= {True, $Failed, ./TCPComm}
```